

Control and application stage	Notes
UNSTAINED CONTROL Include: Unstained sample with no antibody ideally <u>in every experiment</u>	Control for background autofluorescence used to determine positivity - Consider auto fluorescence profile of the sample when choosing flurochromes. - Try using a different excitation laser or detector if autofluorescence is high .
COMPENSATION CONTROL Include: Single stained tubes for each color in the panel ideally <u>in every experiment</u> (attention to tandem lot and uncoupling)	Compensate spillover in multicolor panels where fluorescence overlap between flurochromes is expected (use antibody capture beads or cells) - Use same fluorophore or dye than the sample. In tandems same lot must be used. - Controls must follow the same protocol than the sample (fixation and permeabilization buffers) - Controls must have a positive and negative population with identical autofluorescence - Positive population must be as bright or brighter than in the sample - Collect enough number of events to be statistically significant

Staining Control	
Control purpose and application stage	Notes
VIABILITY CONTROL <i>Staining Control (No Specific Binding): Dead Cells Control</i> Include: DNA based dyes or amine reactive dyes ideally <u>in every experiment</u>	Remove staining artifacts and no specific binding (NSB) from highly autofluorescent dead cells - Use DNA cell impermeant dyes for nucleic acid stain in cells with compromised membranes - Use anime reactive dyes (fixable dyes) for fixed cells or to increase detection options.
REAGENT CONTROL <i>Staining Control (No Specific Binding, NSB): Electrostatic Interactions</i> Include: antibody titrations <u>at characterization stage</u> or new lot introduced	Reduce NSB due to electrostatic interactions - Too high concentrations cause more non-specific binding - Especially important in intracellular antigen staining
FC BLOCKING CONTROL <i>Staining Control (No Specific Binding): Surface Fc receptors</i> Include: FC blocking reagents in protocol prior to staining <u>at characterization stage/ in every experiment</u>	Avoid increased background and false positives due to Ab binding though cell surface Fc receptors - Saturate Fc receptor with specific Abs (FcBlock CD16/32) or preincubate with the serum of the host of the primary antibody - Use of commercially available Fab fragment (REAffinity) antibodies
ISOTYPE CONTROL <i>Staining Control (No Specific Binding)</i> Include: cells labeled with abs with irrelevant specificity (same isotype and flurochromes) <u>at characterization stage</u>	Assess the effectiveness of the Fc blocking protocol The use of isotype controls for gating purposes is not recommended (antibody and isotypes may show different affinities and different F:P ratios)
SECONDARY ANTIBODY CONTROL <i>Staining Control (No Specific Binding)</i> Include: Cells stained only with secondary Ab <u>at characterization stage or new lot introduced</u>	Evaluate nonspecific binding of the secondary antibody Only for samples with indirect staining with secondary antibody or streptavidin

Analysis Control	
Control and application stage	Notes
FMO CONTROLS <i>Analysis Control</i> Include: cells stained with all flurochromes minus one for each color in the panel <u>in every experiment</u>	Detect fluorescence spread in a particular detector introduced by the rest of the flurochromes in the panel. Aids in gating placement. - Especially important in the case of dimly expressed antigens or when antigen expression is unknown. - Also in detectors where fluorescence spillover from bright flurochromes is expected.