


Aurora CS sorter Unidad de Citometría de Flujo CNIC

Ultraviolet laser

355 nm




16 channels for detection 365-829 nm

Detector names

UV1–UV16

Violet laser

405 nm




16 channels for detection 420-829 nm

V1–V16

Blue laser

488 nm




14 channels for detection 498-829 nm

B1–B14

Yellowgreen laser

561 nm




10 channels for detection 567-829 nm

YG1–YG10

Red laser

640 nm



8 channels for detection 652-829 nm

R1–R8