

FACSCanto II
 Unidad de Citometría CNIC

DETECTOR	LASER	PARAMETROS									DICROICO	BANDPASS	RANGO
B	405	V 450/50	Alexa Fluor 350	Alexa Fluor 430	Pacific Blue	CFP	DAPI	bechst 332	FL7			450/50	425_475
A	405	V 510/50	AmCyan	FL8							505LP	525/50	500_550
E	488	B 530/30	Alexa Fluor 488	FITC	GFP	YFP	YFP	FL1			505LP	530/30	515_545
D	488	B 585/42	PE	PI	FL2						556LP	586/42	564_606
B	488	B 670LP	PE-Cy5	PerCP-Cy5.5	PerCP	PI (fl3)	FL3				655LP	670LP	670_735
A	488	B 780/60	PE-Cy7	PE-H7	FL4						735LP	780/60	750_810
C	633	R 660/20	APC	Alexa Fluor647	Cy5-5	CyChrome	FL5					660/20	650_670
A	633	R 780/60	APC-Cy7	APC-H7	Alexa Fluor 750	FL6					735LP	780/60	750-810